

Product Information

Potassium Chloride Solution (0.075 M)

Cat. No. KCL-B (100 ml)

General Information

Potassium Chloride (0,075M) is a sterile-filtered, hypotonic solution commonly used to enlarge cells for an adequate spreading of metaphase chromosome for karyotyping. The solution is prepared in cell culture grade water.

Applications

- Karyotyping

Product Specifications

Appearance	Clear liquid
Storage and shelf life	Store at +2°C to +8°C
CAS No. (solid)	7447-40-7

Precautions and Disclaimer

This product is for research use only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).