

## Technical specifications of TripleA 1000 µl Tips

Applicable products: WB 5064

### Characteristics:

Volume	1000 µl
Low retention	Yes
Sterile	Yes
Filter	No
Packaging	Rack with hinge
Graduations	At 100, 500, and 1000 µl
Tip opening	Tapered for full sample release
Color	Colorless
Compatibility	i.a. Biohit, Brand, Eppendorf, Thermo/Finn, Gilson, Rainin, Socorex

For further information on product sterility and material composition please refer to the appendix

- COA Westburg tripleA tips
- Material composition sheet Westburg-plastics

**The dimension and design of the 1000 µl tip including the position of the filter. Dimensions are in millimeters.**

