

Technical specifications of TripleA 200 µl Tips

Applicable products: WB 5040 - WB 5140 – WB 5145 - WB 5540 - WB 5640

Characteristics:

Volume	200 µl
Low retention	Yes
Sterile	Yes (WB 5145)
Filter	Yes (WB 5040), No (WB 5140, WB 5540, WB 5640)
Packaging	Rack with hinge (WB 5040, WB 5140, WB 5145), Bulk (WB 5540), Eco refill (WB 5640)
Graduations	At 10, 50 and 100 µl
Tip opening	Tapered for full sample release
Color	Colorless
Compatibility	i.a. Biohit, Brand, Eppendorf, Thermo/Finn, Gilson, Rainin, Socorex

For further information on product sterility and material composition please refer to the appendix

- COA Westburg tripleA tips
- Material composition sheet Westburg-plastics

The dimension and design of the 200 µl tip including the position of the filter. Dimensions are in millimeters.

