

Technical specifications TripleA 1250 µl Tips

Applicable products: WB 5074 - WB 5174 – WB 5178 – WB 5570 - WB 5670

Characteristics:

Volume	1250 µl
Low retention	Yes
Sterile	Yes (WB 5178)
Filter	Yes (WB 5074), Nee (WB 5174, WB 5570, WB 5670)
Packaging	Rack with hinge (WB 5074, WB 5174), Bulk (WB 5570), Eco refill (WB 5670)
Graduations	At 100, 500, and 1000 µl
Tip opening	Tapered for full sample release
Color	Colorless
Compatibility	i.a. Biohit, Brand, Eppendorf, Thermo/Finn, Gilson, Rainin, Socorex

For further information on product sterility and material composition please refer to the appendix

- COA Westburg tripleA tips
- Material composition sheet Westburg-plastics

The dimension and design of the 1250 µl tip including the position of the filter. Dimensions are in millimeters.

