



**PurePrep Mini** - Compact Automated Nucleic Acid Extraction System  
(1-16 samples simultaneously, UV lamp sterilization)

**REDUCE HANDS-ON TIME AND  
GUARANTEE HIGH YIELDS**

The PurePrep Mini reduces hands-on time and increases productivity in lab research by automating the purification of nucleic acids, doing away with tedious manual extraction. With its small footprint saving precious workspace and mobile app allowing you to set programs remotely, the system is a must-have for small research environments. By utilizing 6-well tube strips, the PurePrep Mini even allows you to run single samples without wasting plastic consumables.

The innuPREP extraction kits are optimized for usage with the PurePrep extraction systems, allowing extraction from a wide range of sample types.

PurePrep Mini      845-PP-0016000

Available Kits for HMW and UHMW DNA Extraction

innuPREP SE HMW DNA Kit -PP Mini	845-PSS-4216016 /96	845-PSP-4216016/96
innuPREP SE UHMW DNA Blood & Cells Kit - PP Mini	845-PSS-4016016 /96	845-PSP-4016016/96

Available Kits for Magnetic Bead Extraction

innuPREP Genomic DNA Kit - PP Mini	845-PS-0020016	845-PS-0020096
innuPREP Anipath DNA/RNA Kit - PP Mini	845-PS-0030016	845-PS-0030096
innuPREP FFPE DNA Kit - PP Mini	845-PS-0040016	845-PS-0040096
innuPREP Forensic DNA Kit - PP Mini	845-PS-0050016	845-PS-0050096
innuPREP RNA Kit - PP Mini	845-PS-0060016	845-PS-0060096
innuPREP Blood DNA Kit - PP Mini	845-PS-0070016	845-PS-0070096
innuPREP Plate Set - PP Mini	845-PS-1000016	845-PS-1000096
innuPREP Stripe Set - PP Mini	845-PS-2000016	845-PS-2000096

**For more information contact our team: [info.innu@ist-ag.com](mailto:info.innu@ist-ag.com)**