

## Description

- Serological pipettes 1ml white box in single bag sterile

## Part number

- 923000

Pieces by sales : 1000

---

## Description

- Volume: 1 ml
- Material: Polystyrene
- Product manufactured in a facility that is registered to the current version of ISO 9001 or ISO 13485
- Sterile (irradiation)
- Non-Pyrogenic
- DNase/RNase free - Endotoxin free
- BSE/TSE
- Volumetric Accuracy - Serological pipets are accurate to  $\leq 2\%$  at full volume in compliance with ASTM E93
- Packaging: 100 individually sterile wrapped pipettes per bag

