

Description

- Serological pipettes 2ml white box in single bag sterile

Part number

- 923002

Pieces by sales : 800

Description

- Volume: 2 ml
- Material: Polystyrene
- Product manufactured in a facility that is registered to the current version of ISO 9001 or ISO 13485
- Sterile (irradiation)
- Non-Pyrogenic
- DNase/RNase free - Endotoxin free
- BSE/TSE
- Volumetric Accuracy - Serological pipets are accurate to $\leq 2\%$ at full volume in compliance with ASTM E93
- Packaging: 100 pipettes individually packed per bag - 25pcs/zip bag

