

Description

- Serological pipettes 50ml white box in single bag sterile

Part number

- 923008

Pieces by sales : 100



Description

- Volume: 50 ml
- Material: Polystyrene
- Product manufactured in a facility that is registered to the current version of ISO 9001 or ISO 13485
- Sterile (irradiation)
- Non-Pyrogenic - DNase/RNase free - Endotoxin free
- BSE/TSE
- Volumetric Accuracy - Serological pipets are accurate to $\leq 2\%$ at full volume in compliance with ASTM E93
- Packaging: individual packed - 25pcs/zip bag

