

Description

- Serological pipettes 10ml white box in single bag sterile

Part number

- 923010

Pieces by sales : 200

Description

- Volume: 10 ml
- Material: Polystyrene
- Product manufactured in a facility that is registered to the current version of ISO 9001 or ISO 13485
- Sterile (irradiation)
- Non-Pyrogenic - DNase/RNase free - Endotoxin free
- BSE/TSE
- Volumetric Accuracy - Serological pipets are accurate to $\leq 2\%$ at full volume in compliance with ASTM E93
- Packing : Individual packed - 50 pcs/zip bag

