

Product Information

Chicken Serum, sterile-filtered

Cat. No. CHI-1A (500 ml); CHI-1B (100 ml)

General Information

Chicken Serum is a nutrient-rich fluid derived from the blood of healthy chickens without any red blood cells or other clotting components. It is used as a supplement for cell culture media and as blocking agent in immunoassays or as an antibody diluent.

Product Specifications

Appearance	Amber liquid
Storage and shelf life	Store at $\leq -15^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened, store at $+4^{\circ}\text{C}$ and use within 4-6 weeks.
Shipping conditions	Frozen (Dry ice)
Thawing	Overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$. Swirl gently to homogenize.

Precautions and Disclaimer

The product is for research use only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).