

Product Information

CHO|ONE E

Expression Medium, w/o Insulin, w/o L-Glutamine, with Pluronic™

Cat. No. CHOE-500ML (500 ml), CHOE-1000ML (1000 ml)

General Information

CHO|ONE E allows the cultivation of recombinant CHO cells for fed-batch cultures in spinner flasks or large-scale bioreactors. It can be used for effective cultivation of recombinant CHO cells, achieving high yield of the active protein.

CHO|ONE E already leads to increased cell growth and protein yield but when combined with CHO|ONE Feed 1 and Feed 2 the full potential of the Media System can be deployed.

Product Specifications

Appearance	Clear yellow orange solution
Storage and shelf life	Store at +2°C to +8°C protected from light. Once opened store at 4°C and use within 6-8 weeks.
Shipping conditions	Ambient

Formulation

This formulation is our proprietary composition and has no counterparts either in its composition, or in its action.

Precautions and Disclaimer

This product is for research use and further manufacturing only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).